

# CERTIFIED®

A handy and portable spotting kit that can be used on any size spotting job.

# SPOT & STAIN

## Rescue Kit



### Kit Contains

- Quart size bottles with pour tops
- Fresh Spot (W-1)
- Coffee / Tea (W-3)
- Digester
- Dr. Spot (Aerosol)
- White Towel
- Tamping Brush
- Spotting Chart

C361-004

### SPOTTING CHART

Problem	Fresh Spot W-1	Coffee-Tea W-3	Digester	Dr. Spot
Adhesives				
Alcohol Drinks				
Animal Glues				
Asphalt				
Auto Grease				
Beer				
Blood				
Beverage Stains				
Fresh Beverage				
Candy				
Chewing Gum				
Clay/Dirt				
Chocolate				
Coffee/Tea				
Coffee/Tea w/cream				
Cola				
Cosmetics				
Crayon				
Cream				
Feces				
Fruit Juices				
Greasy Foods				
Ice Cream				
Inks (washable)				
Liquor				
Lipstick				
Mildew				
Milk				
Motor Oils				
Mustard				
Paints (latex)				
Paint				
Saliva				
Perfumes				
Shoe Polish				
Soft Drinks				
Sugar				
Syrup				
Urine				
Tar				
Vomit				

All Certified's spotting agents are ready-to-use. All of the above problems can be removed with one or more spotters.

### WHEN USING A SPOTTING KIT

- 1- Identify your spot, and reference on the chart above.
- 2- Pre-test products in inconspicuous location.
- 3- Apply spotter by blotting and tamping.  
Do not rub!
- 4- Always work from edge to center of spot.
- 5- Repeat on area until problem is gone.
- 6- Rinse with water. Do not over wet.

UnoClean

WWW.UNOCLEAN.COM

TOLL FREE: 1-888-226-2724

