

OB-1600

# ObboMed® Medical Air Mattress

Suitable for Stage 1-4 Pressure Ulcer

## ObboMed® UtilityAir™ Mattress

6" High Alternating • 7" High Side Bolster • Low Air Loss • 15° Self-Lateral Rotation

**6" DEEP  
1 LAYER MATTRESS**

**15° SELF-ROTATION**

**EASY TO OPERATE  
CONTROL PANEL**

**7"H SIDE BOLSTER**



The ObboMed® difference—

# Comfort, Safety & Ease




User friendly control panel basic functions:

- On/Off (Standby)
- Lock
- Max Inflate (Auto Firm)
- Static Mode
- Weight and Pressure Display
- Low Air Loss up to 15 LPM
- Audio and Visual Warning Signals (if power outage and/or when air pressure is too low)
- Fowler (Seating) is designed to increase support in the seat/ torso section of the mattress. When the head of the bed is raised, the fowler is intended to prevent any potential risk of patients from "bottoming out."
- Alternating Cycle Time: 10/20/30/40 Min (system defaults to "Alternating" mode automatically after Max inflate or 1-hour in Static Mode)
- Rotation Circle: 10/20/30/40 Min

Features of ObboMed's multi-layered mattress:

- Designed to minimize friction and shearing of the skin
- Vapor permeable to facilitate moisture removal into the airflow away from the mattress
- Perspiration vacates through the mattress cover surface and/or the patient's bed boundary
- Antimicrobial, Liquid resistant, Fireproof
- Latex, lead, and mercury-free
- Below the cover layer, vapor flows out with the Low Air Loss feature to lower humidity under the patient while simultaneously transferring away moisture and heat
- Air cells are composed of a combination of thermoplastic polyurethane (TPU) and nylon and heat-melt airtight with ultrasonic high-frequency, and are liquid proof and machine washable



Patented UtilityAir™ Features\*:

**15° Self-lateral Rotation:** This feature has specially designed continuously alternating air cells creating a gentle, wave-like motion that provides a subtle rotation. This therapy aids in redistributing the user's weight to alleviate pressure points while stimulating blood flow. The subtle motion also helps improve sleep.

\* For details on these patented features, visit [obboMed.com/patents](http://obboMed.com/patents)

Additional UtilityAir™ Features\*:

Side Bolster

Sinking feature lowers to 4" deep to keep patient centralized and prevent falling off bed during rotation

OB-1600 MATTRESS DIMENSION AND WEIGHT SPECIFICATIONS

1-layer Air Cell, 6"H Mattress, 7"H Side Bolster

**Mattress Dimension:** 80"L x 36"W x 7"H  
**Mattress Weight:** 12 lbs  
**Pump:** 30 LPM  
**Weight Capacity:** 450 lbs

**Control Unit Weight:** 10 lbs  
**Unit Package Weight:** 24.7 lbs  
**Unit Package Dimension:** 18.5"L x 16.5"W x 21.3"H

Mor-Medical International in Partnership with ObboMed Medical Air Mattress

Limited Warranty:

1 year for Top Cover, 2 years for Control Unit and Mattress

©2024 ObboMed Group LLC. All rights reserved. Duplication, republication, or modification in whole or partially is prohibited without prior written permission from ObboMed. All trademarks are owned by or licensed to ObboMed Group LLC or its subsidiaries unless otherwise noted.